

Hyland Laboratories
4501 Colorado Blvd.
Los Angeles 39, Calif.

Please send Tissue Culture catalog to

Name _____

Organization
or Firm _____

Street _____

City _____ Zone _____ State _____



**SEND FOR HYLAND'S
LATEST LISTING OF**

TISSUE CULTURE COMPONENTS

Of special interest is unique *Newborn Agamma Calf Serum*, which provides an excellent protein source for cell propagation and is recommended for detection, propagation and study of many viruses. This specially processed serum, from which gamma-globulin has been completely removed by fractionation technics, provides an unusually high content of alpha- and beta-globulins. In virus studies, it has shown no inhibition of Types I, II and III polio viruses. This serum is available in both liquid and dried form in a variety of practical sizes.

Other bovine specialties include: *Newborn Calf Serum* (liquid or dried), which is derived from 1- to 4-day-old calves and, because of high alpha-globulin content, is more stimulatory to cell growth than serum from more mature animals; non-toxic *Fetal Calf Serum* (liquid); *Bovine Amniotic Fluid* (liquid or dried); *Bovine Embryo Extract* (dried); *Bovine Embryo Extract, Ultrafiltrate* (liquid); *Bovine Serum* (liquid or dried); *Bovine Serum, Ultrafiltrate* (liquid).

Chicken Serum (liquid or dried) is available in large pools for polio testing. Our line also embraces other serums and serous fluids, ultrafiltrates, balanced salt solutions and synthetic media. We welcome your inquiries about special formulations or preparations of particular interest to you.

HYLAND LABORATORIES



4501 Colorado Blvd., Los Angeles 39, Calif.



Microtiter[®]

**8 TIMES AS MANY SEROLOGICAL AND
OTHER DILUTIONS ACCOMPLISHED WITH
THE SAME AMOUNT OF REAGENTS FORMERLY USED**

Our microtitrator kit (MICROTITER)* utilizing individually calibrated spiral loops and pipette droppers permits fast and accurate titrations.

Each kit furnished with two types of plates to permit complement fixation, hemagglutination, hemagglutination inhibition, and metabolic inhibition tests.

You can use eight loops simultaneously . . . and use the kit case as a wash rack for the plates.

**MADE IN AMERICA BY AMERICAN LABOR
AND FROM AMERICAN RAW MATERIALS**

Immediate delivery and prompt service on replacement parts.

* PATENT PENDING

write for free illustrated literature

Cooke Engineering Company

735 North St. Asaph Street, Alexandria, Virginia



● an entirely new concept—simple as 1-2-3

Complete
Tissue
Cultures

at
7.4 pH

ready to use
in seconds



1 CulturSTAT Reagents are supplied in two 100 ml bottles, with a sterile, disposable Vacutainer needle. Bottle 1 contains desiccated perishable reagents of the culture media. Bottle 2 contains a balanced salt solution. CulturSTAT Reagents have long shelf life at normal room temperatures.



2 To restore the desiccated material, insert the short end of the Vacutainer needle through the seal of Bottle 2 and the long end through the seal of Bottle 1. Vacuum automatically transfers the salt solution aseptically in seconds.



3 Withdraw the needle and the complete, sterile tissue culture, at 7.4 pH \pm 0.2, is ready for use in seconds.

CulturSTAT

This product, originated by Cappel Laboratories, eliminates entirely the time and effort required to obtain and pipette individual ingredients of a tissue culture medium and to adjust to the proper pH. The possibility of pipetting errors is eliminated and uniformity is assured by standardization.

Ask your dealer for CulturSTAT Tissue Culture Reagents today.

Cappel, CulturSTAT and Vacutainer are Trademarks



Cappel Laboratories

ANIMAL BLOOD PRODUCTS • TISSUE CULTURE REAGENTS • DIAGNOSTIC REAGENTS

West Chester, Pennsylvania

DIVISION OF B-D LABORATORIES, INC.

For Reliability in
Tissue Culture Reagents

depend on
Cappel

STERILE ANIMAL SERA

Plasma and Embryonic Extracts

CHEMICALLY DEFINED REAGENTS

(Stock and 10X Stock)

Maintenance, Balanced Salt and Special Solutions


TISSUE CELL CULTURES

Serially Propagated and Primary Cultures

CulturSTAT A new, time-saving concept in prepared tissue culture reagents

Prepared with rigid adherence to original formulae, and standardized, Cappel Tissue Culture Reagents provide properties which assure stability and full control of laboratory techniques. Their increasing use by leading laboratories is evidence of product quality, superior value and service.

Specify, through your distributor, drop shipment of fresh Cappel Tissue Culture Reagents from regular production schedules



Cappel Laboratories

ANIMAL BLOOD PRODUCTS • TISSUE CULTURE REAGENTS • DIAGNOSTIC REAGENTS

West Chester, Pennsylvania

DIVISION OF B-D LABORATORIES, INC.

An important new tool
in Bacteriology . . .

Skerman's A GUIDE TO THE IDENTIFICATION OF THE GENERA OF BACTERIA—*With Methods and Digests of Generic Characteristics*

"The Key to the Genera of Bacteria has been compiled with the aim of placing in the hands of research workers, teachers and students, a volume in which general directives for the identification of bacteria are supported by a complete list of the techniques needed for the purpose. In addition, a Digest of data published for the various species in the seventh edition of *Bergey's Manual of Determinative Bacteriology* and in several original papers has been included. The publication of the Key and Methods in the one book should encourage a more general application of common procedures for the description of bacteria."—*Preface*

By V. B. D. SKERMAN, *Reader in Bacteriology and Head of the Department of Bacteriology, University of Queensland, Brisbane, Australia*
1959 • 567 pp., 32 pls. • Paper cover • \$5.50



THE WILLIAMS & WILKINS COMPANY

Baltimore 2, Maryland

THE JOURNAL OF IMMUNOLOGY

Information for Subscribers

Published monthly, beginning January, 1929. (From 1916–1925, and in 1928 Volumes I–X, inclusive, and Vol. XV published bi-monthly.)

The subscription price is \$18.00 per year in U.S.A.; \$18.50 in Canada; \$19.50 elsewhere. Single copies, when available, are \$2.25 in U.S.A.; \$2.50 elsewhere. The subscriptions of members of the American Association of Immunologists are included in the membership dues and are handled by the Association's Treasurer, Dr. Byron H. Waksman, Massachusetts General Hospital, Boston 14, Mass. All other subscriptions should be sent to the publishers.

Change of address: Publisher must be notified 60 days in advance. Journals undeliverable because of incorrect address will be destroyed. Duplicates can be obtained (if available) from the publisher at the regular price of single issues.


New subscriptions and renewals are entered to begin with the first volume of the year when subscription is received. Should any issue of the current volume be out of print at the time the subscription order is received, the pro-rata value of such numbers will be refunded to the subscriber.

Subscriptions should be renewed promptly: To avoid a break in your series, subscriptions should be renewed promptly. The publishers cannot guarantee to supply back issues on belated renewals.

Correspondence concerning all business matters should be addressed to

THE WILLIAMS & WILKINS COMPANY

Publishers of Scientific Books and Periodicals
BALTIMORE 2, MARYLAND, U. S. A.



PROBABLY THE WATER IS ALL RIGHT!

A cool drink from the common cup. A small chance. Usually no one suffers for it. Usually. It's human nature to take chances. Except at Nutritional Biochemicals. There, human nature gives way to perfectionism. Because lives depend on the absolute purity of N.B.Co.'s biochemicals. So does successful research. In ordering one of N.B.Co.'s 2600 biochemicals, you express confidence. Confidence in a company whose only business is preparing pure biochemicals. A company whose world-wide volume brings you pure biochemical

at low prices. Send for our Free catalogue today. Or call us at M0ntrose 2-0214, Cleveland, Ohio.

NUTRITIONAL BIOCHEMICALS CORPORATION
21010 Miles Avenue • Cleveland 28, Ohio
24-Hour Delivery in the U.S.A. • Slightly Longer Anywhere Else

Send for our free October, 1961 Catalog containing more than 2600 items. Fill out coupon and mail today for your copy. J1

Name _____
Organization _____
Address _____
City _____
State _____ Zone _____

