

Prophylactic Effectiveness of Live and Killed Tularemia Vaccines. I. Production of Vaccine and Evaluation in the White Mouse and Guinea Pig. HENRY T. EIGELSBACH AND CORA M. DOWNS.....	415
Mouse Antibody. I. Characterization and Properties of Antibody in Mouse Peritoneal Fluid. R. L. ANACKER AND J. MUNOZ.....	426
Biological Activity of Aggregated $\gamma$ -Globulin. IV. Mechanism of Complement Fixation. TERUKO ISHIZAKA, KIMISHIGE ISHIZAKA AND TIBOR BORSOS.....	433
The Allergens in Ragweed Pollen. Physicochemical and Immunologic Studies on a Purified Fraction of the Dialysate of Ragweed Pollen. M. RICHTER, B. W. GRIFFITHS AND B. ROSE.....	439
Bacterial Adherence: A Method for Detecting Antibody Production by Single Cells. O. MÄKELÄ AND G. J. V. NOSSAL.....	447
Study of Antibody-Producing Capacity of Single Cells by Bacterial Adherence and Immobilization. O. MÄKELÄ AND G. J. V. NOSSAL.....	457
The Inactivation of Rinderpest Virus by $\beta$ -Propiolactone and Its Effect on Homologous Complement-Fixing and Neutralizing Antibodies. S. S. STONE AND P. D. DELAY..	464
Interaction of Rickettsiae and Phagocytic Host Cells. II. Chemotactic Action of Typhus Rickettsiae on Human Polymorphonuclear Leukocytes <i>in Vitro</i> . CHARLES L. WISSEMAN, JR.....	468
Antigenicity, Solubility and Chemical Stability of Renin Complexes. ERWIN HAAS AND HARRY GOLDBLATT.....	472
Antibodies to Denatured Desoxyribonucleic Acid in a Lupus Erythematosus Serum. DAVID STOLLAR AND LAWRENCE LEVINE.....	477
Endotoxin Susceptibility and Delayed Hypersensitivity in Experimental Histoplasmosis. EDITH D. BOX AND NORMAN T. BRIGGS.....	485
Role of the Adrenals in Stress Induced Resistance to Anaphylactic Shock. PERRY E. TREADWELL AND A. F. RASMUSSEN, JR.....	492



## Production experience guarantees **RELIABILITY**

Your laboratory and research needs demand the quality and dependability which Colorado Serum Co. assures its customers. Prompt service. All correspondence and inquiries answered immediately.

- serums
- bloods
- ultrafiltrates
- complement
- globulins
- fluorescent materials
- diagnostic reagents
- tissue culture reagents

We maintain a variety of our own laboratory animals under the finest conditions.

The newest equipment, coupled with over 35 years of production experience and devotion to product quality guarantees you can order with confidence.

Write for this FREE CATALOG NOW!  
No salesman will call.



**COLORADO SERUM CO.**  
**LABORATORIES**

Laboratory and General Office

4950 YORK STREET • DENVER 16, COLORADO • MAin 3-5373



PEAK OF QUALITY

# PPLO

## *ISOLATION and CULTIVATION*

THE AVAILABILITY of media for the isolation and cultivation of *Pleuropneumonia-Like Organisms* (PPLO) has increased the interest in these infectious agents. PPLO are readily grown from clinic specimens or from pureculture on Bacto-PPLO Agar and Bacto-PPLO Broths enriched with Bacto-PPLO Serum Fraction, or Bacto-Ascitic Fluid.

**BACTO-PPLO AGAR** is an especially prepared solid medium containing the essential factors from meat infusion required by the PPLO. This medium, enriched with Bacto-PPLO Serum Fraction, will support abundant growth of all PPLO strains.

**BACTO-PPLO BROTH w/o CV** enriched with Bacto-Ascitic Fluid or Bacto-PPLO Serum Fraction is a liquid medium for enrichment and cultivation of PPLO organisms. Selective media for PPLO may be prepared from Bacto-PPLO Broth w/o CV by the addition of thallium acetate, potassium tellurite, penicillin or other bacteriostatic agents as desired.

**BACTO-PPLO BROTH with CV** enriched with Bacto-Ascitic Fluid is a selective medium for the isolation and cultivation of PPLO from specimens contaminated with commensal organisms. Contains Crystal Violet for inhibition of Gram-positives. The selectivity of the medium may be further increased to inhibit Gram-negatives by the addition of Bacto-Chapman Tellurite Solution. This medium will support growth of all PPLO strains and effectively suppress the growth of most contaminating bacterial flora.

**BACTO-PPLO SERUM FRACTION** contains the growth accessory factors required by all strains of the PPLO and is used as an enrichment in Bacto-PPLO Agar and Bacto-PPLO Broth w/o CV.

**BACTO-ASCITIC FLUID** is a sterile, bile-free ascitic fluid specially prepared for enriching culture media used in the cultivation of the PPLO group and other fastidious microorganisms.

**BACTO-CHAPMAN TELLURITE SOLUTION** to be added to Bacto-PPLO Broth w/CV or Bacto-PPLO Broth w/o CV for the inhibition of Gram-negative microorganisms.

THE DIFCO MANUAL, NINTH EDITION,  
*including descriptions of these media and their use,*  
*is available on request.*

**DIFCO LABORATORIES**  
DETROIT 1, MICHIGAN