



Spot the difference.



Ready to move from one platform to another with complete confidence? Invitrogen's BioSource™ Luminex® Assays are the only ones benchmarked to ELISA and calibrated to NIBSC standards or published intracellular cell models for phosphoprotein activity. We also test standards before and after lyophilization, so you can count on consistent performance. With more than 150 kits to choose from and expert field support ready to work by your side, Invitrogen's Luminex® Instrument and Assays can help deepen your understanding of extracellular and intracellular proteins. See how you can cross over at www.invitrogen.com/luminexassay.

 **invitrogen™**

Welcome to a more colorful worldSM

Dazzling performance with BD FACSAria II

Transform your expectations of multicolor capability

Achieve a new level of productivity with the BD FACSAria™ II cell sorter. Maintain optimal efficiency, sensitivity, and resolution, even at high-speed sort settings. The system makes multicolor more accessible with its fixed-alignment cuvette and new setup and tracking software.

The BD FACSAria II opens up a whole new world of sorting possibilities. As always, multicolor-optimized BD reagents, consulting services, and application support are standard. Call your BD representative or visit bdbiosciences.com/aria2 to learn more.



Helping all people
live healthy lives

Class I (1) laser product.
This product is for Research Use Only. Not for use in diagnostic or therapeutic procedures.
BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2008 BD
23-10597-01

BD Biosciences
2350 Qume Drive
San Jose, CA 95131
bdbiosciences.com